ENCISION recommends placing this document in the Instructions for Use/Care Section of your AEM Monitor Operator/Service Manual.

Device Description

End Point Monitor Remote Display

The End Point Monitor Remote Display is an optional accessory for use only with the AEM® Monitor to provide a remote display of the bipolar current indication.

Principle of Operation

The End Point Monitor Remote Display consists of a 30 segment bar graph display. The numbers of segments illuminated in the bar graph display corresponds to the magnitude of the bipolar current as indicated on the scale underneath.

NOTE

- This product is not intended to be located in the sterile field. It is not supplied sterile.

Use with Bipolar Surgery

In bipolar electrosurgery, the surgical instrument includes both electrodes. A patient return electrode need not be used. Current flows between the two electrodes and through the tissue contacted by the instrument, heating the tissue. In bipolar electrosurgery, close control is needed to ensure the correct degree of heating.

End Point Monitoring of the bipolar instrument will assist the surgeon in confirming the end point of bipolar desiccation. This information is displayed on the End Point Monitor Remote Display as an illuminated visual bar graph.

Indications/Intended Use

For use with the AEM Monitor when remote monitoring of bipolar current is desired.

Contraindications

None known.

Instructions For Use

Prior to Use

Thoroughly read these instructions and the instructions in the AEM Monitor Operator/Service Manual.

Inspect the instrument for proper assembly and function.

AEM System Setup/Assembly Instructions

See laminated Setup Sheet (00701) when using the ES9005 series AEM Cord Adapter and (02678) when using the ES9015 Universal Adapter.

1. Connect the End Point Monitor Remote Display into the receptacle marked Bipolar End Point Monitor Remote Display located on the rear panel of the AEM Monitor.
2. Mount the End Point Monitor Remote Display near the TV monitor or at another location in the view of the operating staff.

NOTE

- When the End Point Monitor Remote Display is connected, the AEM Monitor bar graph display does not function.
- Do not allow fluids to enter the End Point Monitor Remote Display.
- The End Point Monitor Remote Display is not designed to withstand dropping onto a hard floor. Route the display cable to avoid the likelihood of tripping.

CAUTION

- Good operating room practice suggests that connections of accessories to electrosurgical generators be made only while the generator is Off or on Standby.
- See electrosurgical generator manual and AEM Monitor Operator/Service Manual for precautions concerning the general application of Bipolar Surgical Equipment.

Cleaning Instructions

Follow the procedures approved by your institution or use a validated infection control procedure.

Use a mild cleaning solution of disinfectant and a damp cloth to thoroughly wipe all surfaces and the connecting cable.

WARNING

- Electric shock hazard. Always turn off and unplug the AEM Monitor and accessories, including the End Point Monitor Remote Display, before cleaning.

CAUTION

- Do not clean the End Point Monitor Remote Display with abrasive cleaning or disinfectant compounds, solvents or other materials that could scratch the panels or damage the unit.
- Do not steam sterilize the End Point Monitor Remote Display.

Technical Specifications

Duty Cycle:

Continuous

Bipolar Current Indicators Function (RF Amperes):

Accuracy: 10% of full scale

Range: 0.1 to 1.1 Ampere

Maximum current without damage: 3.0 Amperes

Dimensions:

1.89 in. (7.34 cm) wide x 5.26 in. (13.26 cm) long x

0.79 in. (2.00 cm) high

Weight with Cord:

1.0 lb. (0.45 kg)

Environmental Characteristics:

Operating Temperature: 50 to 122° F (10 to 50° C)

Storage Temperature: -13 to 158° F (-25 to 70° C)

Atmospheric Pressure (Operating/Storage): 50 to 110 kPa

Express Warranty

ENCISION hereby warrants to Buyer that products purchased hereunder shall be free from defects in materials and workmanship under normal use and service, as specified in this Instructions for Use/Care, for the period of six months from date of purchase. Any evidence of unauthorized repair, modification or sterilization of this product will void this warranty.

See AEM Monitor Operator/Service Manual for details of Limitations, Disclaimer, and Exclusions.

Return of Used Product

If for any reason this product must be returned to ENCISION, a returned goods authorization is required prior to shipping. Appropriate return instructions may be obtained from ENCISION.

Product

ENCISION reserves the right to amend, modify or to change any product, to introduce new products, to withdraw products and otherwise vary product specifications at any time without notice.

ENCISION® and AEM® are registered trademarks of ENCISION Inc.

	Federal (USA) law restricts this device to sale by or on the order of a physician.
	Consult Instructions for Use
	Latex Free

Made in USA

 **Manufactured by**
ENCISION Inc.
6797 Winchester Circle
Boulder, CO 80301 USA
Ph: 303-444-2600
Fax: 303-444-2693

EC	REP
----	-----

Authorized Representative
(according to MDD93/42/EEC)
MDSS GmbH
Schiffgraben 41
30175 Hannover, Germany

Printed in USA
© Copyright 2010 Encision Inc.

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ENCISION® en AEM® zijn gedeponeerde handelsmerken van ENCISION Inc.